

About the Purify O3

Purify O3: The Simple, Effective Way To Clean Your CPAP/BiPAP Mask & Equipment



Purify O3 uses ozone (activated oxygen) technology to sanitize your equipment, so you can rest easy knowing your CPAP/BiPAP equipment is germ and bacteria free! Ozone disinfection is one of the most effective methods to sanitize areas of PAP equipment that cleansing wipes and UV light cannot reach.

Why Sanitize With Ozone Technology?

Natural: Activated oxygen, known as ozone (O₃), is a naturally formed gas that is waterless and chemical free. Ozone is widely used as a disinfecting agent in a number of industries world-wide, including medical facilities.

Safe & Effective: Purify O3 provides hospital grade sanitation to quickly, permanently kills 99.99% of viruses and bacteria.

Affordable: Purify O3 includes all you need to sanitize your device in one affordable kit. No extra pre-wash solutions, filters or replacement bags required. HSA/FSA eligible.

Simple & Easy To Use

One-Touch Controls: One button push-to-start initiates the disinfection process. The unit will shut off automatically when the cycle is complete.



Quick Clean: Disinfection cycle lasts only 35 minutes, with an additional 90 minute recommended rest to allow residual ozone to dissipate.

Maintenance Free: Liquid-free cleaning. No additional pre-wash solutions or filters required.

Ultra Quiet: Ultra quiet sanitizing, with audible alerts when cycle is complete.

Travel Easy: Lightweight and portable design saves space and travels easily alongside your PAP machine. The long-lasting battery allows multiple cleanings from one charge. Dual purpose sanitizing/travel case included in kit.



Don't Take a Chance With Your Health!

FDA Registered!