



*The Safe, Effective Method to Sanitize G5 Machines and Applicators*

**KILLS 99.99%  
OF GERMS  
& BACTERIA**



## Specifications

**Size** - 4.5" x 2.75" x 2" (11.4 x 7 x 5 cm)

**Weight** - 0.5 lbs (0.227 kg)

**Cleaning Technology** - Ozone (O3)

**Power Source** - 4W, Rechargeable Lithium-Ion Battery

**Charging Source** - USB (cable and AC charging block included)

**Includes** - Purify O3 unit, dual purpose sanitizing and storage bag, USB cable & AC charging block, instruction manual

## Simple To Use:

- Place G5 Applicator or Neocussor into Sanitizing Bag. For larger G5 machines, place Purify O3 and G5 machine into a sealed plastic bin.
- Closed container sanitizing process ensures safe and effective equipment sanitization at a higher efficacy level than other cleaning methods.
- Push the easy one-touch control to start. Automatically shuts off after 35 minute cleaning cycle is complete.
- Recharges quickly using wall outlet or USB port.
- No maintenance, filters, or extra cleaning solutions required. Nothing additional to purchase ... ever!
- May also be used to clean CPAP/BiPAP & nebulizer masks, breast pump accessories, and other medical devices.



**FI Sales, LLC®**

St. Charles, IL 60174 USA  
(630) 444-1380  
Fax: (630) 444-1381

Visit FI Sales online at  
[www.fisalesllc.com](http://www.fisalesllc.com) or  
[info@fisalesllc.com](mailto:info@fisalesllc.com)