



# Flimm-Fighter®

## Technical Specifications



Motor -	1/35 horsepower
Output Speed -	20-30 cycles per second
Power Source -	110v or 220v (other voltages available)
Speed Control -	Solid State Electronic
Current Leakage -	under 75 MicroAmps
Power Cord -	10' hospital grade
Drive Cable -	5'
Housing Material -	ABS Cicolac
Configurations -	Table-top mounted
Unit Dimensions -	10" H x 4" W x 4" D
Unit Weight -	8 Lbs.
Shipping Dimensions -	22" x 21" x 17 <sup>3/8</sup> "
Shipping Weight -	11 Lbs. (varies based on applicator package)
Applicator Packages -	Includes AP250, AP251*, AP253, AP234 (Self Application Package, SAPKG, sold separately) *For the Pediatric Flimm-Fighter, AP252 is substituted for AP251



**GENERAL PHYSIOTHERAPY, INC.**  
13222 LAKEFRONT DRIVE - EARTH CITY (ST. LOUIS), MISSOURI 63045-1504 - USA  
PHONE: (314) 291-1442 TOLL FREE: (800) 237-1832 FAX: (314) 291-1485  
[www.g5.com](http://www.g5.com)

